



Multiplex LC-MS/MS lysosphingolipids analysis in plasma for the screening of sphingolipidoses and Niemann-Pick type C disease


Magali PETTAZZONI

Service de Biochimie et Biologie Moléculaire Grand Est
Unité des Maladies Héréditaires du Métabolisme
Centre de Biologie Est
Hospices Civils de Lyon
FRANCE



Monique Piraud, Cécile Pagan, David Cheillan, Christine Vianey-Saban, Roseline Froissart.

Sphingolipidoses




The measurement of the primary accumulated sphingolipids in plasma is poorly informative

Diagnosis based on enzymatic testing

Filipin staining test

Sphingolipids/LysoSphingolipids in plasma



From Sun and Zhang,
Advances in Gaucher Disease: Basic and Clinical Perspectives 260 pages,
August 2013. Future Medicine Ltd.


- Lysosphingolipids (LysoSL) = deacylated form of the sphingolipids
- Role in pathophysiology ?
- With the use of MS/MS Emerging biomarkers **in plasma** for screening of sphingolipidoses and Niemann-Pick type C

Tandem Mass Spectrometry

**Very sensitive,
very specific technique**

Allows **the simultaneous measurement** of numerous compounds **all together**, in complex mixtures, according to **specific fragmentation of each molecule**

« **TRANSITION** » = pair of (precursor ion)/(product ion)



LC-MS/MS LysoSLs measurement

LysoGlobotriaosylceramide (LysoGb₃) Fabry

Aerts et al. Proc Natl Acad Sci USA, 2008 (not MS/MS)
Boutin et al. Clin Chim Acta, 2012

LysoGalactosyl-ceramide (LysoGalCer) Krabbe

Chuang et al. Clin Chim Acta 2013

LysoGlucoyl-ceramide (LysoGlcCer) Gaucher

Dekker et al. Blood 2011

LysoSphingomyelin (LysoSM) Niemann-Pick A/B

Chuang et al. Mol Genet Metab 2014

LysoSM analogue 509 (LysoSM509) Niemann-Pick A/B and C

Giese et al. Orphanet J Rare Dis 2015

~~**LysoSulfatide (LysoSulf) Metachromatic leukodystrophy**~~

Mirzaian et al. Blood cell Mol Dis 2015

DE GRUYTER Clin Chem Lab Med 2016; 000

Open Access

Giulia Polo, Alessandro P. Burlina, Thilini B. Koламunnage, Michele Zampieri, Carlo Dionisi-Vici, Pietro Strisciuglio, Martina Zaninotto, Mario Piabani and Alberto B. Burlina*

Diagnosis of sphingolipidoses: a new simultaneous measurement of lysosphingolipids by LC-MS/MS

Polo et al. Clin Chem Lab Med 2016

Constantinoble available on ScienceDirect

Clinica Chimica Acta

Journal homepage: www.elsevier.com/locate/clinchem

Simultaneous quantitation of sphingoid bases by UPLC-ESI-MS/MS with identical ¹³C-encoded internal standards

M. Mirzaian¹, P. Wissar², M.J. Ferraz³, A.R.A. Marques⁴, P. Gaspar⁴, S.V. Gussarsen⁴, K. Kytidou⁴, J.D.C. Codée⁴, C. van der Marf⁵, H.S. Overkleeff¹, J.M. Aerts^{1,6*}

Mirzaian et al. Clin Chim Acta 2017

RESEARCH ARTICLE

LC-MS/MS multiplex analysis of lysosphingolipids in plasma and amniotic

+ LysoGM1 ganglioside (LysoGM1)

+ LysoGM2 ganglioside (LysoGM2)

Magali Pettazoni^{1*}, Roseline Froissart^{1,2}, Cécile Pagan¹, Marie T. Vanier^{3,4}, Séverine Ruet¹, Philippe Latour⁵, Nathalie Guffon⁶, Alain Foulhoux⁶, Dominique P. Germain⁷, Thierry Lavade⁸, Christine Vianey-Saban^{1,9}, Monique Piraud^{1,9}, David Cheillan^{1,9}

Pettazoni M et al. PLoS One. 2017

Method

Extraction/Purification

200 µL of EDTA plasma + IS + MeOH
SPE (Oasis MCX, 30 mg, 60 µm, Waters Corp)

Liquid chromatography

C8 column (Uptisphere. Interchim©)
Gradient of elution (Phase A: H₂O 0.2% FA, Phase B ACN 0.2% FA)

MS/MS– Multiple Reaction Monitoring Mode (MRM)




API 4500 Q-Trap
(AbSciex©)

Analyte	STANDARDS	Internal Standard **	MODE	TRANSITION		
Lyso Gb3	Lyso Gb ₃	Glycinated lysoGb ₃	Positive	786.5	>	282.3
Lyso HexCer	Lyso GlcCer	Lyso GlcCer d5	Positive	462.3	>	282.3
Lyso SM	Lyso SM d18:1	Lyso SM d17 :1	Positive	465.3	>	184
Lyso SM 509	Lyso SM d18:1	Lyso SM d17 :1	Positive	509.3	>	184
Lyso GM1	Lyso GM1	S1P d17 :1	Negative	1278.6	>	290.1
Lyso GM2	/	S1P d17 :1	Negative	1116.6	>	290.1*

* No standard; calculated from mass/structure

** Commercially available

Positive mode



LysoGLUCOSYLceramide and lysoGALACTOSYLceramide: not separated
Measured as lysohexosylceramide « LysoHexCer »


Quantitative validation
Except LysoGM1 and LysoGM2

LysoSL multiplex measurement:
Routine experience in Lyon
for the screening of sphingolipidoses and NPC

Fabry disease (X-linked)

Mutations on *GLA* gene

Classical forms / Variant forms with specific variant mutations on *GLA* gene



MALES DIAGNOSIS

α Gal A activity (always deficient)

Fabry females

Heterozygote females

- *Normal* α Gal A: 30% cases
- *Normal* (U)Gb₃ :
18% cases (classical forms)
most of cases (variant forms)

FEMALES DIAGNOSIS


- α Gal A activity
- (U) Gb₃
- *GLA* gene

Fabry disease: Plasma LysoGb₃

Results
LYON,
France

		Males Classical form	Males Variant form*	Females Classical form	Females Variant form**	Controls
	N = 77	n = 15	n = 10	n = 44	N = 8	N = 228
LysoGb ₃ (nmol/L)	Mean	91	5.8	5.7	1.4	< 0.6 (99 th perc.)
	Range	30 - 173	1.8 - 11.1	1.1-31.2	0.3-4	0-1.1

Normal in
4/8 cases





* Variant *GLA*
mutations in males

p.Ile232Thr
p.Met296Val
p.Arg301Gln
p.Arg363His
p.Asn215Ser (n=2)
p.Asn215Ser
p.Ile198Thr
p.Phe113Leu (n=3)

** Variant *GLA*
mutations in
females

p.Arg363His (n=3)
p.Asn215Ser
p.Phe113Leu (n=4)


Fabry disease: Correlation (U)Gb₃ / (P) LysoGb₃



Gaucher disease: LysoHexCer

Results LYON, France

		GAUCHER	Saposin C deficiency *	KRABBE (infantile)	KRABBE (juvenile / adult)	Controls
		n = 27	n = 1	n = 8	n = 4	n = 228
LysoHexCer (nmol/L)	Mean	160	76	13	1.2	< 3.3 (99 th perc.)
	Range	5 - 427	-	9 - 22	0.9 - 1.4	0.1-3.5




* Treated by SRT

- Highly elevated in GD, large range of values, can be very moderately elevated in some cases
- Elevated in case of saposin C deficiency
- In GD: Correlation with chitotriosidase and CCL18 (= indirect and non specific biomarkers)

Gaucher disease: LysoHexCer/*GBA* Genotype correlation


Biological variability of LysoHexCer levels +++ for a same genotype

Results of the French cohort, untreated Gaucher patients



Lower LysoHexCer values in patients homozygous for N370S mutation of *GBA* gene

Gaucher disease: LysoHexCer in Amniotic Fluid (AF)




In case of non-immune hydrops foetalis (NIHF)
New powerful biomarker in AF of fetal GD

		Controls	Gaucher AF	NIHF Controls AF
		n = 15	n = 5	n = 64
LysoHexCer (nmol/L)	Mean	< 0.7	460 ***	< 0.7
	Range		68 - 996	

Krabbe disease: LysoHexCer

Results
LYON,
France

		GAUCHER	KRABBE (infantile)	KRABBE (juvenile / adult)	Controls
		n = 27	n = 8	n = 4	n = 228
LysoHexCer (nmol/L)	Mean	160	13	1.2	< 3.3 (99 th perc.)
	Range	5 - 427	9 - 22	0.9 - 1.4	0.1-3.5



Moderate increase in infantile KD
In normal range in juvenile and adult forms
 Saposin A deficiency ? To be evaluated

Niemann-Pick diseases

Niemann-Pick type A/B

Sphingomyelinase deficiency

- Tedious enzyme activity measurement

Niemann-Pick type C

Lysosomal cholesterol trafficking and lipid storage disorder

- Filipin staining test on cultured fibroblasts
- Genetic testing *NPC1* or *NPC2* gene

Chitotriosidase

Moderate but inconstant elevation

6 % of patients are homozygous for 24 bp duplication in *CHIT1* gene and 30 – 40 % heterozygotes

MS/MS

Oxysterols plasma (2010)

Cholestane-3 α ,5 α ,6 α -triol (and 7-ketocholesterol)

LysoSM (2014)

Molecular Genetics and Metabolism 111 (2014) 392–9

Contents lists available at ScienceDirect

Molecular Genetics and Metabolism

Journal homepage: www.elsevier.com/locate/ymgme

Short Communication

Lyso-sphingomyelin is elevated in dried blood spots of Niemann–Pick B patients

Wei-Lien Chuang^a, Joshua Pacheco^a, Samantha Cooper^a, Margaret M. McGovern^a, Gerald E. Cox^a, Joan Seutzer^a, X. Kate Zhang^{a,*}

^a Center for Pediatric and Family Genomics, The Miriam Hospital, Providence, RI, 02912-0178, USA

^b Department of Pediatrics, Brigham Young University, Utah Valley University, Provo, UT, 84602, USA



LysoSM « 509 » = LSM + 44 (2015)

Giese et al. *Orphanet Journal of Rare Diseases* (2015) 10:78
DOI: 10.1186/s13023-015-0174-1

ORF ORPHANET JOURNAL OF RARE DISEASES

RESEARCH Open Access

A novel, highly sensitive and specific biomarker for Niemann–Pick type C1 disease


Anne-Karin Giese^{1†}, Hermann Mascher^{2†}, Ulrike Grüner³, Sabrina Bucher⁴, Guido Ramp⁴, Jan Lukas¹, Danielle Vučković⁵, Naci Al Eise⁵, Mario Corina-Borja⁶, Forbes D Porter⁷, Frances M Platt⁸ and Anett Rolfs^{1,4*}



NPA/B and NPC: LysoSM and LysoSM509 in plasma

Results LYON, France

		NPA/B	NPC* < 10 y	NPC* > 10 y	Controls
		n = 24	n = 26	N = 29	n = 228
LysoSM (nmol/L)	Mean	16.5	0.7	1.0	< 1.9 (99 th perc.)
	Range	2.4 - 69.6	0.2 - 2.1	0.2 - 3.5	0.1-2.0
LysoSM509 (MoM)	Mean	214	149	63	< 4.3 (99 th perc.)
	Range	98-515	34 - 390	7 - 144	0-9.1




* Some patients are treated


GM1 & GM2 gangliosidoses: LysoGM1/GM2 in plasma

IS: d17 Sphingosine-1P

LysoGM1



LysoGM2



In spite of poor recovery results (lack of specific IS)

Not detected in controls

Abnormal presence of lysoGM1
and lysoGM2

in GM1 gangliosidosis

in Sandhoff diseases

in Tay-Sachs diseases

Good specificity

(5/6 cases)

(7/9 cases)

(3/4 cases)

Lack of sensitivity
for mild and adult forms

Conclusion

Our LysoSLs multiplex assay: 6 lysoSLs including LysoGM1 and LysoGM2

- **Efficient and rapid biochemical screening tool**
 - ✓ Small plasma sample volume (pediatrics +++)
 - ✓ Better screening of FD males with variant form, and females
 - ✓ More sensitive screening of NPC than oxysterols
 - ✓ Differential screening of NPC, NPA/B and GD in the same run
 - ✓ Specific screening of fetal GD
 - ✓ Useful for the monitoring of patients (FD, GD)
- **Positive results to be confirmed by enzymatic measurement, and molecular studies**
- **Adults and mild phenotype: possibly normal screening**
- **Other analogues of LysoSL** (mass variation on the sphingosine moiety): could be interesting biomarkers to be evaluated ?
(C Auray-Blais, Sherbrooke, Canada)

Acknowledgments



Unité des Maladies Héréditaires du Métabolisme
Centre de Biologie et de Pathologie Est

Monique Piraud

Roseline Froissart

Cécile Pagan

David Cheillan

Séverine Ruet

Christine Vianey-Saban

Marie-T. Vanier

Philippe Latour

Alain Fouilhoux

Nathalie Guffon

Dominique P Germain

Thierry Levade

Michèle Trauchessec

Florine Coeurdray

Virginie Gaucher



Hospices Civils de Lyon

magali.pettazzoni@chu-lyon.fr